

Reve

Restore Your Sleep

- CPAP Re-engineered
- Personalized Comfort
- Supportive Connected Care

Reve is a **next-gen CPAP device** integrating advanced comfort features, connectivity, and safety—all while being whisper-quiet without any noise-reducing foam, thanks to patented honeycomb technology.



**Auto
CPAP**

FEATURES

Packed for Comfort, Compliance & Clinical Excellence

- **Therapy Modes**
 - 1. Fixed CPAP Mode
 - 2. Auto PAP Mode
- **Advanced Event Detection**
 - 1. Obstructive Apneas & Hypopneas
 - 2. Central Apneas & Hypopneas
- **Smart Comfort & Sleep Initiation**
 - 1. Timed Ramp & Smart Ramp (iRamp): Adjusts ramp timing based on real sleep detection
 - 2. Comfort Control+ (CC+)
 - 3. Intelligent Start: Device starts when the patient breathes into the mask

SPECIFICATIONS

Pressure Range	4 to 20 cmH ₂ O (0.5 cmH ₂ O increments)
Sound Level	Only 24.5 dB at 10 cmH ₂ O
Max Turbine Flow	190 l/min at 19 cmH ₂ O
Dimensions	176 x 144 x 104 mm (with side cover) 234 x 144 x 104 mm (with humidifier)
Weight	1.0 kg (with side cover) 1.3 kg (with humidifier)
Data Collection	SD card, USB, Bluetooth
Oximetry	Optional sensor integration

 **MADE IN FRANCE**

