

Reve

Auto CPAP

Restore Your **Sleep**

CPAP Re-engineered

Personalized Comfort

Supportive Connected Care



Restore Your Sleep
with **Reve**

MANUFACTURED BY

Sanrai International

Head office: 1 Heritage Circle Suite 6, Endicott, NY 13760, USA

Phone: +1 607 821-0153 Email - info@sanrai.com

www.sanrai.com

MARKETED BY

Sanrai Med India

Head office: Office Tower B, Elante, Office No. 506 Plot No.178-178A

Industrial Area Phase-1 Chandigarh- 160002

Phone: +91 73 4110 3371 Email - info@sanraimed.com

www.shop.sanraimed.com

**MADE IN FRANCE**

ISO 0459 CE



CPAP RE-ENGINEERED



Foamless

Uses innovative patented honeycomb technology for your safety



24.6 dB, Whisper Soft

Sleep uninterrupted with the quietest CPAP available



190 L/min flowrate

Effortlessly maintains open airways for all body sizes



Circuit Select

Calibrates your circuit resistance to deliver more effective therapy



Mask Fit & Go

Measures the leak of each mask and seal to calibrate pressure for more efficient therapy



Detects Central and Obstructive Apneas


Uses cardiac oscillations to detect both obstructive and central apneas and hypopneas



Honeycomb technology is a novel way to separate the air flowing through **Reve to reduce noise without using foam.**

Designed

for better Sleep



PERSONALIZED COMFORT



Auto CPAP

Automatically adjusts pressure to maintain open airways for effective therapy



Intelligent Start and Stop

Reve automatically starts or stops when you put on or remove your mask



Intelligent Ramp

Gradually increases pressure after detecting that you have fallen asleep



Comfort Control+

3 comfort levels gently reducing pressure making it easier when you exhale



Heated Humidification

10 settings to customize humidification levels to warm your air and prevent a dry mouth



Heated Tubing

Automatic temperature control to prevent condensation in your tubing while asleep



With advanced comfort settings, getting accustomed to Auto CPAP therapy has never been easier. With the **personalized settings of Reve, being compliant is a breeze.**

Comfort

redefined. Rest restored

Reve

SUPPORTIVE CONNECTED CARE



User Friendly

Easy to follow on-screen instructions with an intuitive touchscreen



Patient Support Team

We are here to guide you and optimize your therapy



Connected Therapy

Cellular, Bluetooth, and Wifi connectivity providing real-time therapy insights and telemonitoring



Intuitive and Customizable dashboard

Multiple parameters for clinicians and homecare providers to offer ongoing treatment support



Expert Settings

Includes reducing pressure and setting limits for specialized cases



Access Lite

Phone app to support your treatment journey and monitor progress



With connected care, providers can easily follow along a patient's nightly therapy, **allowing for timely interventions to improve compliance.**

Dedicated

Patient Support Team